Wiptam®

Cobalt-based dental wire for clasps, retainers and friction-pins

User instructions
 Gebrauchsanleitung
 Gebruiksaanwijzing
 Instructions d'utilisation
 Istruzioni per l'uso
 Instrucciones para el uso

Wiptam® is a part of the Wisil® system

Cobalt-based dental wire for clasps, retainers and friction-fit pins. EN

WIPTAM is a colour-fast and mouth-resistant wire whose specific characteristics provide considerable advantages when fabricating clasps, retainers and friction-fit pins. WIPTAM wires are supplied medium hard (F max. 1,000 - 1,200 MPa). When WIPTAM wires are cold-worked, e.g. during bending they reach a spring hardness. Unlike FeCrNi clasp wires, no difficulties will occur when soldering WIPTAM wire. The hardness that is lost due to the heat applied during soldering is regained during mechanical post-treatment (bending, trimming or polishing). WIPTAM wire can also successfully be used for orthodontic objects that require soldering.

Indications for use:

Cobalt-based dental wire for clasps, retainers and friction-pins.

Contraindications:

WIPTAM wire contains nickel. It should not be used by or for people with a known Nickel allergy, or with an allergy or sensitivity to one of the other components.

WIPTAM wire is not suitable for casting channels.

WIPTAM wire cannot be used for casting.

Warnings

In the mechanical processing of dental alloys, always work with local aspiration plus face or respiratory protection.

Cobalt - vapour is released during processing. May cause cancer.

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May release metal vapours when melting. Do not breathe dusts and vapours. Ensure suitable suction/aeration at the work place and with operational machinery.

Adverse effects:

Allergies to metals contained in the alloy as well as electrochemically conditioned dysesthesia are possible. In individual cases systemic side effects of metals contained in the alloy have been reported.

Interactions:

Avoid occlusal and approximal contact of different alloy types.

Package:

The package contains eight straightened 25-cm wire rods.

Using WIPTAM:

To ensure the stability of the clasp basis, always use wire that is thicker than the clasp to be made. Before soldering, bend the clasp with smooth-tipped pliers. Sharp-edged tools or pliers with rough beaks may create grooves, which are potential weak spots. White gold solder or WITEX CoCr solder produce the best soldering results. After soldering, adapt the clasp optimally along the contours of the element. Cut the wire to the desired length and taper the end of the clasp. Bending, trimming and polishing WIPTAM wire changes its material structure, converting it from medium hard to springy hard. The greater the force exerted to bend and fit the wire, the greater the increase in wire hardness. If the clasp has become too hard at one particular spot due to intensive mechanical treatment, this specific spot can be heated to make it flexible again.

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