FINO

Keramiktiegel

41380 / 41381 / 41390 / 41391 / 41400 /  
 41408 / 41415 / 41430 / 41433 / 41571 /  
 41572

Keramik-Schmelzmulden  
 41369 / 41370 / 41388 / 41419 / 41428

Deutsch

### Gebrauchsinformation

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#### User Information

##### 1. Fields of application

For melting and casting of dental alloys. With a specifically developed soft inside lining for a long service-life. The additional use of melting powder and crucible protection is not required. Suitable for all common casting units with induction heating or open flame casting.

##### 2. Advantages

- Optimally protected crucible walls
- High number of casts per crucible
- Less residual melt
- Less contaminations of the melt
- Increased quality of the cast

##### 3. Safety precautions

Do not terminate cast process after the melt is already liquid – there is a risk that the crucible fuses.

Before first use the ceramic crucibles and casting pits are heated in the furnace; possible transport damages become thus visible.

In case of heavy cinder deposits the crucibles should be exchanged.

##### 4. Processing instructions

4.1 Information about the inside lining  
Please do not remove the special lining. The soft lining does not react with the melts.

Please do not glaze crucibles and pits.

4.2 Removal of melt residue, cinder  
Do not strain the crucible unnecessarily neither mechanically, thermally nor chemically. To avoid fine cracks in the structure the crucible should not be beaten. Never use a pointed object to scratch cinder or residual melt out of the crucible. This could damage the lining and the chemical attack on the crucible wall will be facilitated. Ceramic is brittle breaks relatively easy.

Residual melt and cinder can often easily removed by using compressed air when the crucible is still hot.

4.3 Information about heating and cooling  
Crucible and melt should not be heated if not necessary.

Extended melting time boosts formation of cinder and increases corrosion. The service-life of the crucible is extended if it is not quenched but left to cool on a sand layer or something similar. Pre-heat the crucible before casting.

The best is to add it to the moulds in the pre-heating furnace. When the crucible is well heated the melt will flow out better.

##### 4.4 Placement of cast alloy

By placing the alloy cones it is important to understand that the metal expands during heating. Wedged metal pieces could therefore break the crucible open.

##### 4.5 Use of flux

Do not glaze crucible with flux. Flux softens every ceramic crucible material. All FINO Ceramic Crucibles have already a special lining. In case the manufacturer of the alloy makes the use of a flux compulsory it must be used extremely sparingly.

##### 4.6 Information for the use of graphite crucibles

Avoid long heating times of the graphite crucibles. The casting unit is heated to the required temperature and then the graphite crucible is inserted.

If installed the lid of the casting unit is now closed.

After the graphite crucible is red-hot the alloy is placed in the crucible according to the directives of the manufacturer, melted and then cast as quickly as possible.

After the casting process is finished leave the switched-off casting unit cooling with inserted graphite crucible.

If the cooling process must be fast we recommend cooling in a sand bath.

Therefore place the crucible in a bowl and cover with sand.

Strictly do not cool in water because this will cause initially invisible cracks in the graphite body.

Never leave the graphite crucible in a heated casting unit without reason because any additional time in a heated casting unit boosts burn-off.

##### 5. Storage

Store in a dry place.

##### 6. Technical data

No further data.

##### 7. Delivery forms

**FINO Ceramic Crucibles**

10 pcs. 41380

For casting units Heraeus Kulzer. Types CL-G, CL-G 97 and Heracast RC S.

10 pcs. 41381

For casting units Heraeus Kulzer. types CL-G 77, CL-G 94 and Heracast RC L.

5 pcs. 41390

**FINO Ceramic Crucibles special shape for Fornax**

(For casting unit Bego, type Fornax.  
special shape: The bevelled shape of the crucible accounts for a faster flow of the melt and less residues after the cast).

5 pcs. 41391

For casting unit Bego, type Fornax.

5 pcs. 41400  
For casting units Heraeus Kulzer, types CL-IG, CL-IM, CL-I 95, Heracast EC and iQ.

5 pcs. 41408  
For casting unit Dentsply DeguDent, type Degutron.

5 pcs. 41415  
For casting unit Linn, types Minitherm, HFS-3 Standard.

5 pcs. 41430  
For casting unit Schütz, type Microtronic.

5 pcs. 41433  
For casting unit Manfredi, type Multihertz (old version).

1 pc 41571  
5 pcs. 41572

For casting units Galloni, types Modular 3N, Modular 3S and all Fusus.

#### **FINO Ceramic Casting Pits**

5 pcs. 41369  
For casting units Dentsply DeguDent, types Motorcast compact and Motorcast TS 1.

1 pcs. 41370  
For casting units Dentsply DeguDent, types Motorcast compact and Motorcast TS 1.

5 pcs. 41388  
For casting units Bego, types Fandor, Castor and Pollux. (deep).

5 pcs. 41419  
For casting unit Kerr, type Universal.

5 pcs. 41428  
For casting units Reitel, types Retocast T and Retocast TE.

#### **8. Guarantee**

Our technical recommendations of application are based on our own experiences and tests and should only be regarded as guidelines.

It rests with the skills and experience of the user to verify that the products supplied by us are suitable for the intended procedures. Our products are undergoing a continuous further development. We reserve the right of changes in construction and composition. It is understood that we guarantee the impeccable quality of our products.

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