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**GUARANTEED TOP QUALITY
THROUGHOUT**

English

User Information

FINO FLOWCAST Ceramic Crucibles

41591 / 41592

By purchasing this product you have made an excellent choice.

FINO quality products convince through consistently high quality for economical prices.

1. FIELDS OF APPLICATION

For the melting and casting of dental alloys.

2. ADVANTAGES

- ◆ Crack resistant
- ◆ Extremely high resistance against temperature changes
- ◆ Long service life

3. PHYSICAL DATA

Pre-heating temperature before first usage:
800 °C

4. PROCESSING / INSTRUCTIONS

Do not overstrain the crucible neither mechanically nor thermal nor chemically. To avoid structural cracks the crucible should not be dusted. Never rasp the crucible with sharp instruments to remove slag or residues of the charge. This will damage the coating and facilitates chemical damage to the wall of the crucible. Crucible and charge should not be heated-up unnecessarily. Extended

melting time promotes the formation of slag and boosts corrosion.

The service life of the crucible is extended, if the crucible is not cooled down rapidly but can cool slowly on a layer of sand. The crucible should be pre-heated before casting. It is recommended to place it together with the moulds in the pre-heating furnace. A uniformly heated-up crucible ensures a smooth ejection of the charge. If the crucible shows a heavy deposit of slag it should be discarded in time.

When the crucible is loaded it should be remembered that the alloy undergoes a strong expansion during heating-up. Wedged pieces of alloy can crack the crucible. Do not glaze the crucible with flux. Flux will macerate any ceramic material crucibles are made of. All FINO ceramic crucibles come with a special coating. The special coating is glazed prior to the first use at 800 °C / 5 min.

5. STORAGE

Store in a dry place.

6. SAFETY PRECAUTIONS

Do not interrupt the casting process when the charge is already liquid – it is an imminent hazard that the crucible is melted through.

7. AVAILABILITY

FINO Ceramic Crucibles 5 pcs.

8. GUARANTEE

Our technical recommendations of application are based on our own experiences and tests and should only be regarded as guidelines. It rests with the skills and experience of the user to verify that the products supplied by us are suitable for the intended procedures. Our products are undergoing a continuous further development. We reserve the right of changes in construction and composition. It is understood that we guarantee the impeccable quality of our products.

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