REF 42009/42009P

CE 0483 ISO 9333

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

CE 0483 ISO 9333

ENG Instructions for use

Please read carefully before use
Keep for referral purposes

1. Intended purpose

CoCr - solder for soldering dental restorations.

2. Indication

For soldering cobalt chrome partial denture alloys. Specifically adapted to the requirements of the FINOLLOY partial denture alloy.

- FINOLLOY H
- FINOLLOY XH

Not for ceramic veneering.

2.1 Contraindication

Do not use if incompatibilities or allergies to constituents of the alloy have been confirmed.

3. Fields of application

Only for processing in the dental laboratory according to the intended use.

4. Safety precautions

Danger

H317 May cause an allergic skin reaction.

H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.

H413 May cause long lasting harmful effects to aquatic life.

 Date of manufacture

12.2 Application

Please read carefully before use
Keep for referral purposes

Rx only

Only to be used by qualified personnel!

Do not reuse

Biohazard

12.3 Storage

Temperature limit

Keep away from direct sunlight

Protect against moisture

Do not use damaged packaging

top - store upright

Caution fragile

4.1 Prevention

P273 Avoid release to the environment.

P285 In case of inadequate ventilation wear respiratory protection.

P363 Wash contaminated clothing before reuse.

4.2 Response

P302+P352 IF ON SKIN: Wash with plenty of water.

P304+P341 IF INHALED: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing.

P342+P311 If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician.

4.3 Side effects

No side effects known. Individual or allergic reactions to constituents of the alloy or electrochemically-induced discomfort are possible in rare cases. In this case, do not use the product.

4.4 Interactions

In case of occlusal or interproximal contact of different alloys electro-chemically related irritations may be possible in rare individual cases.

4.5 Precautions

There is no information available on the safety or efficacy for treatment of children and pregnant or nursing women. FINOLLOY may cause interference with the evaluation of MRI examinations.

Wear dust mask!

Metal dust is a health hazard. For blasting and processing use extraction and dust mask type FFP3-EN149:2001.

5. Processing instructions

5.1 Preparation

Grind areas for soldering.

For soldering repairs, widen the defective area in the shape of a funnel.

Create a solder block from solder investment material (e.g. FINOTHERM) and pre-dry in the pre-heating furnace.

If a solder block is not required, the object should be clamped in a solder holder.

5.2 Soldering

When soldering in the solder block, keep the gap to be soldered approx. 0.2 mm to max. 0.5 mm wide. The soldering points must be exposed so that they are easy to access with the flame from all sides.

Apply a fluxing agent to the solder joint, break off an appropriate piece of FINOLLOY CoCr Solder and coat with flux.

Then apply the solder and heat the solder joint with a circling flame.

The solder can also be tack welded.

Coat with fluxing agent before soldering.

As soon as the object to be soldered glows red, point the flame at the solder until it fills the solder joint.

The outline of the piece to be soldered is still visible.

6. Storage

P402 Store in a dry place.

7. Physical data

Processing temperature: 1180 °C; Melting interval: 1125 - 1195 °C.

7.1 Composition in %:

Co	Cr	Si	Mo	Fe	B
61,0	28,5	4,2	3,1	1,3	1,5

Additional elements < 1%: C.

Free from nickel, cadmium, beryllium and lead.

8. Delivery forms

FINOLLOY CoCr solder sticks	4,5 g	42009
Trial pack	1 g	42009P

9. Accessories

FINO DEOXYD flux	80 g	42060
------------------	------	-------

10. Guarantee

Our technical recommendations of application are based on our own experiences and tests and should only be regarded as guidelines. It rests with the skills and experience of the user to verify that the products supplied by us are suitable for the intended procedures. Our products are undergoing a continuous further development. We reserve the right of changes in construction and composition. It is understood that we guarantee the impeccable quality of our products.

11. Disposal

P501 Dispose of contents/container in accordance with local/regional/national/ international regulations.

11.1 EWC waste code product
170407 CONSTRUCTION AND DEMOLITION WASTES (INCLUDING EXCAVATED SOIL FROM CONTAMINATED SITES); metals (including their alloys); mixed metals.

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