Original FINO  
Germany

# FINOFORNO S FINOFORNO L

Vorwärmofen  
Pre-heating Furnace  
four de préchauffage  
Horno de precalentamiento

43186

43187

FINO · der feine Unterschied · the fine difference · la fine différence · la fina diferencia · la differenza fine

Gebrauchsinformation · User Information · Mode d'emploi  
Información sobre el uso · Informazioni d'uso  
Gebruiksaanwijzing · Instrukcja stosowania

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


### 10.5 Registrierungsnummer

WEEE-Reg.-Nr. DE 75542610

Version 1/2010. Read the user manual before you start operating the unit!

# TABLE OF CONTENTS

<b>1. General</b>	<b>3</b>
• Information about this user manual • Explanation of symbols • Limitation of liability	
• Copyright protection • Spare parts • Warranty information • Customer service	
• Declaration of conformity	
<b>2. Safety</b>	<b>4</b>
• Responsibility of the user • Staff requirements • Proper use • Personal protective equipment	
• Specific hazards • Measures in case of hazards and accidents • Labelling	
<b>3. Technical data</b>	<b>7</b>
• General information • Connection values • Operation conditions • Type label	
<b>4. Design and function</b>	<b>8</b>
• Overview • Assembly description • Connections • Danger areas • Control elements	
• Delivery includes • Options	
<b>5. Transport, packing and storage</b>	<b>13</b>
• Safety information for transport • Symbols on packaging • Transport inspection	
• Packaging • Transport	
<b>6. Installation and initial operation</b>	<b>14</b>
• Installation • Safety • Installation • Connection to mains supply • Initial operation	
<b>7. Operation</b>	<b>16</b>
• Safety • Loading of furnace • Switch on furnace • Programming • Select program	
• Start program • Start program automatically • Temperatur stepping • Switch off	
<b>8. Malfunctions</b>	<b>22</b>
• Safety • Malfunction indications • Table of malfunctions	
• System information for malfunction diagnosis	
<b>9. Maintenance</b>	<b>24</b>
• Safety • Maintenance schedule • Maintenance work	
<b>10. Disassembling and disposal</b>	<b>25</b>
• Safety • Disassembling • Disposal • Disposal according to EC directive 2002/96/EC	
• Registration number	

## 1. General

### 1.1 Information about this user manual

This user manual enables the safe and efficient operation of the unit.

The manual is part of the unit and must be kept in the direct vicinity accessible by the personnel at all times. The personnel must carefully read this manual and must have understood it before any operation may commence. The basic requirement for a safe operation is the observance of all given safety precautions and operational information in this manual. Furthermore the local accident prevention regulations and the general safety regulations for the field of application of this unit apply.

Pictures in this user manual are meant for the basic understanding and may differ from the actual design of the unit. Apart from this manual the user directions of the integrated components that are attached in the appendix apply.

### 1.2 Explanation of symbols

#### Safety precautions

In this user manual safety precautions are indicated by symbols. The safety precautions commence with signal words, which express the scope of the risk. To avoid accidents, personal injuries and tangible damage unconditionally observe safety precautions.

#### Tips and recommendations

#### NOTE!

Points out helpful tips and recommendations as well as information for an efficient and trouble-free operation.

#### Specific safety precautions

To highlight specific risks the following symbols are used in safety precautions:

#### DANGER! Risk of fatal injury from electrical current!

Indicates perilous situations caused by electric current. If the safety precautions are not observed there is the risk of serious injuries or death.

## DANGER!

Indicates an imminently dangerous situation, which might be fatal or might lead to serious injuries if not avoided.

## CAUTION!

Indicates a possibly dangerous situation, which might cause minor or slight injuries if not avoided.

## CAUTION!

Indicates a possibly dangerous situation, which might cause minor or slight injuries if not avoided.

## CAUTION !

Indicates a possibly dangerous situation, which might lead to tangible damage if not avoided.

## 1.3 Limitation of liability

All information and references in this user manual were compiled taking into account the valid standards and regulations, the state of technology as well as our knowledge and experience of many years. The manufacturer does not assume liability for damages due to:

- Non-observance of the user manual
- Any use not according to specifications
- Operation by non-trained personnel
- Arbitrary alterations
- Technical alterations
- Use of non-certified spare parts

The actual delivery may vary from the here described information and design in case of special models, the application off addition order options or due to most recent technical alterations.

The agreed obligations from the supply contract, our general terms and conditions as well as the supply conditions of the manufacturer and the valid legal regulations at the time of conclusion of contract apply.

We reserve the right to implement technical alterations concerning the improvement of usage requirements and further developments.

## 1.4 Copyright protection

This user manual is protected by copyright and exclusively intended for internal use. Without written authorisation from the manufacturer it is prohibited to surrender this manual to third parties, to make copies of any kind and form – even in parts – as well as use and/or communicate the contents except for internal use. Non compliance obligates indemnity. The right to claim additional indemnities remains reserved.

## 1.5 Spare parts

### **WARNING!** Safety risk due to wrong spare parts!

Faulty or non-genuine spare parts may interfere with the safety and may as well lead to damages, malfunctions or total loss. Therefore:

- Use only genuine spare parts from the manufacturer.

Procure spare parts only from an authorised dealer or directly from the manufacturer.

## 1.6 Warranty information

Warranty information and general terms and conditions can be seen in the internet and can be downloaded (address » last page).

## 1.7 Customer service

For technical advice please contact our customer service (address » last page). Furthermore our staff is always keen to learn about new information and experiences that result from practical application and that may be valuable for the improvement of our products.

## 1.8 Declaration of conformity

The declaration of conformity is included with the product on delivery.

## 2. Safety

This section gives an overview over the most important safety aspects to ensure an optimum protection for the personnel as well as for a safe and trouble-free operation.

The non-observance of the operational information and the safety precautions from this manual may lead to considerable risks.

## 2.1 Responsibility of the user

The unit is employed in the commercial sector. The user of this unit therefore is subject to legal obligations concerning occupational safety. Besides the occupational safety information in this user manual the current safety, accident prevention and environmental protection regulations for the employment of this unit must be observed. Thereby particularly applies:

- The user must familiarise himself with the applicable occupational safety regulations and must additionally establish risks, which are related to specific work conditions at the place of employment of the unit in a risk assessment. This must be implemented in the form of operational directives for the operation of the unit.
- The user must clearly regulate and determine the responsibilities for installation, operation, maintenance and cleaning.
- The user must take care that all personnel that is handling the unit have read and understood the user manual.

Furthermore he must train the personnel and inform about risks in regular intervals.

- The user must provide his personnel with the required protective clothing. Furthermore the user is responsible that the unit is always in a technically perfect condition.

Therefore the following applies:

- The user must take care that the maintenance intervals described in this user manual are adhered to.
- The user must take care that all safety applications for operability and integrity are regularly checked.

## 2.2 Staff requirements

### 2.2.1 Qualifications

#### **WARNING! Risk of injuries due to insufficient qualification!**

Improper use may lead to considerable personal or tangible damage. Hence:

- II operations may only be performed by qualified personnel.

In the user manual the following qualifications are specified for different areas of activity.

#### – Personnel

##### • Qualified laboratory worker / qualified dental technologist

is due to the professional training, knowledge and experience as well as the knowledge of relevant regulations able to perform the assigned work and is also able to self-reliantly identify and avoid possible risk.

##### • Certified electrician

is due to the professional training, knowledge and experience as well as the knowledge of relevant standards and able to perform the assigned work and is also able to self-reliantly identify and avoid possible risk.

## 2.3 Proper use

The unit is exclusively designed and engineered for the here described use:

The pre-heating furnace is exclusively used for the extraction of wax from casting moulds.

### **WARNING!** Danger due to improper use!

Any use of the unit not according to specifications and/or any altered application may lead to dangerous situations. Therefore:

- Use the unit only according to specifications.
- Strictly adhere to all information of this user manual.

The sockets on the back side are exclusively for the connection of:

- a vapour extractor,
- a catalyst
- an extraction hood

Any claims based on damages caused by improper use are excluded. Only the user is liable for all damages caused by improper use.

## 2.4 Personal protective equipment

During work it is required to wear personal protective equipment to minimise health hazards.

- For the respective work always wear the required protective equipment.
- Follow the instructions regarding personal protective equipment that is displayed in the work area.

## Wear for specific work

For the performance of specific work special protective equipment is required. In individual sections of this manual this equipment will be separately referred to. Below this specific protective equipment is discussed:

### Light breathing protection

to protect against harmful dusts.

### Protective gloves

To protect the hands against friction, abrasion, puncture or deeper injuries as well as contact with hot surfaces.

## 2.5 Specific hazards

Please note the here listed safety precautions and warning notes in the further sections of this manual to reduce health hazards and to prevent dangerous situations.

### **DANGER!** Risk of fatal injury from electrical current!

The imminent risk of fatal injury exists if live parts are touched. Damage to insulation or individual parts may be fatally dangerous.

Therefore:

- In case of damaged insulation switch off mains supply immediately and prompt repairs.
- Work on the electric system must only be performed by certified electricians.
- For all work at the electric system switch power off and check that current is off.
- Before any maintenance, cleaning and repair work switch off mains supply and secure against reconnection.
- Do not bridge circuit breakers or put out of operation. If circuit breakers are exchanged keep to the correct Ampere value
- Keep humidity away from live parts. It could cause a short-circuit.

**Harmful vapours****CAUTION! Health risk due to outgassing and vapours during heating operation!**

During heating of waxes and investment materials wax vapours and ammonium can accumulate in the vicinity of the laboratory furnace. Therefore:

- Take care of sufficient ventilation or vapour extraction to outdoors.
- Use vapour extractor DG2 or extraction hood DU1 - DU3/3 with outdoor air exhaust.
- Use catalyst KN1 or KN2.

**Hot surfaces****CAUTION!****Risk of burns due to hot surfaces!**

Contact with hot components may cause burns. Therefore:

- During all work close to hot components strictly wear protective occupational clothing and protective gloves.
- Before any work starts make sure that all components have cooled down to ambient temperature.

**Hot workpieces****WARNING! Risk of burns due to hot workpieces!**

Workpieces reach high temperatures during processing and may cause burns on contact. Therefore:

- Before handling workpieces check if these are hot. If necessary leave to cool.
- Handle hot workpieces only with tongs.

**Mineral fibres**

The heat insulation of the furnace contains parts made from aluminium silicate fibre. This fibre is declared a „carcinogen substance according to category 2 according EU-directive 97/69/EU“ according to EU-classification from 05.12.1997.

In the case of strong and extended dust exposure diseases of the lungs and the pleura such as fibrosis and cancer may be provoked as results from animal experiments show. These results are not confirmed through examinations in humans. Health risks are not to be expected if the recommended use and the applicable threshold values are observed.

Examples for European threshold values:		
Country	Threshold value	Source
Germany	0,5 fibres/ml	TRGS 900
France	0,6 fibres/ml	Circulaire DRT No. 95-4 du 2.01.95
England	2,0 fibres/ml	HSE – EH40 – Maximum Exposure Limit

**2.6 Measures in case of hazards and accidents****– Preventive measures**

- Always be prepared for accidents and fire!
- Keep first aid equipment (first aid kit, blankets, etc.) and fire extinguisher at hand.
- Familiarise personnel with accident reporting, first aid and rescue provisions.
- Clear access roads for rescue vehicles.

**– Measures in case of accidents**

- Activate emergency stop immediately.
- Initiate first aid measures.
- Evacuate persons from the danger area.
- Inform the responsible person on site.
- Call rescue services.
- Clear access roads for rescue vehicles.

**2.7 Labelling**

The following symbols and labels are installed on the outside of the furnace housing:

**WARNING!  
Risk of burns due to illegible symbols!**

With time labels and plates may become soiled or unrecognisable. Therefore:

- Keep all safety and operation labels in a well readable condition.
- Replace damaged plates and labels immediately.

**Warning note "Hot surface"  
(at the FINOFORNO furnace door)**

Hot surface are not always perceivable. Leave to cool respectively wear protective gloves before handling.

**Warning note  
"Electric current"**  
**at the furnace cable duct.**

**3. Technical data****3.1 General information**

Furnace type	FINOFORNO S	FINOFORNO L
External dimensions: (W x H x D) [cm]	40,0 x 48,0 x 40,0	48 x 58 x 52
Chamber dimensions: (W x H x D) [cm]	15 x 10 x 17	20 x 11 x 25
Space for casting moulds:	4 pieces (6x)	6 pieces (9x)
max. temperature:	1100 °C	1150 °C
Heating times:	60 Min. 900 °C 90 Min. 1050 °C	60 Min. 800 °C 120 Min. 1100 °C
Weight:	31,0 kg	70,0 kg


**3.2 Connection values****Electrical**

Furnace type	FINOFORNO S	FINOFORNO L
Voltage:	230 V AC	
Frequency:	50 Hz	
Allowed deviation from set voltage:	± 10 %	
max. Leistungsaufnahme:	1,6 kW	2,3 kW
Consumption / cycle:	Connection to a separate circuit with slow-blow circuit breaker 16 A.	
Fuse protection (by user):	IP 20 (unit is protected against the penetration of foreign bodies but not protected against the penetration of water)	

**3.3 Operation conditions****Environment**

Specification	Value	Unit
Temperature range:	5 - 40	°C
Humidity, maximal:	Up to 31 °C: 80 Up to 40 °C: 50 no condensation	%
Installation:	Indoors in dry rooms	
Height:	max. 2000	m
Degree of pollution:	2	


### 3.4 Type label

Pic . 1: Type label pre-heating furnace

The type label is on the right housing side and contains the following information:


- Manufacturer, address of manufacturer
- Type
- Serial number (S/N)
- Voltage and frequency
- Power consumption
- Disposal information
- CE-certification
- Year and country of manufacture

There is a type label on the optionally available vapour extractor, catalyst and/or additional ventilator respectively that looks similar to the one attached to the pre-heating furnace and that shows information applicable to the option.

## 4. Design and function

### 4.1 Overview

Pic. 2: Sectional view of the pre-heating furnace

- |                   |                                      |
|-------------------|--------------------------------------|
| 1 Controller      | 5 Thermo element (not shown)         |
| 2 Furnace door    | 6 Outlet for exhaust air (not shown) |
| 3 Ceramic insert  | 7 Air circulation blower (option)    |
| 4 Heating chamber | 8 Heat insulation                    |

### 4.2 Assembly description

#### Heating chamber

The heating chamber consists of a high-quality ceramic, which is evenly heated from four sides with heating coils.

In the FINOFORNO S the heating coils are mounted on heating muffles. In the FINOFORNO L they are integrated into the ceramic plates.

The heating coils are covered with chamotte. A ceramic insert protects the heating muffle against contamination.

#### Heat insulation

The heating chamber is from all sides covered with a ceramic insulating layer. In the insulation there is an opening leading to the back side of the unit for the exhaust of vapours.

#### Furnace door

The furnace door has a ceramic door stone for heat insulation. In the FINOFORNO S the door swings open towards the front. In the FINOFORNO L the door slides to the left side. The furnace door is equipped with a safety switch that interrupts the heating current if the door is opened.

#### Controller


The FINOFORNO controller ensures that the desired set-temperature is reached in the shortest possible time. The set-temperature is maintained with an accuracy of 3 °C. A sensor with PtRh-Pt-thermo element inside the heating chamber records the chamber temperature close to the objects to be sintered. A thermo element breaking cut-out prevents the furnace from overheating in case of a sensor defect.

The start time of operation is automatically calculated so that the heating process is completed at the desired time and the sintered objects ready for casting can be removed. Operation parameters and heating programs are saved in a permanent memory and remain saved even in case of a powercut.

#### Furnace housing

The light grey (RAL 7035) furnace housing is made of steel plate and is inside and outside resin-coated.

### 4.3 Connections

Pic. 3: Connection on the back side of the furnace


- 1 Mains connection
- 2 Socket "Gebläse"
- 3 Feed cable thermo element

## 4.4 Danger areas

### Danger area

- Risk of bruising at door mechanic and door crack

**FINOFORNO S**


Pic. 4: Bruising area swing furnace door

**FINOFORNO L**

Pic. 5: Bruising area swivel furnace door

- Risk of burns at the door surface


Pic. 6: Hot surface furnace door

## 4.5 Control elements

### 4.5.1 Furnace door


For pre-heating furnace

**FINOFORNO S**

Pic. 7: Furnace door, Swing mechanism

**FINOFORNO L**


Pic. 8: Furnace door, Swivel mechanism

1 Control element door lever

1 Control element door handle

## 4.5.2 Controller, electronics

Pic. 9: Overview control elements

- |   |                   |   |                  |
|---|-------------------|---|------------------|
| 1 | Operating display | 5 | Arrow keys       |
| 2 | Program level     | 6 | Start/Stop-key   |
| 3 | Program number    | 7 | Casting time key |
| 4 | Main switch       | 8 | Day/time display |

## 4.5.3 Switch and key functions

Pic. 10: Mains switch on		Mains switch lights up if unit is switched on (switch position down)
Pic. 11: Key		Change display mode
Pic. 12: Increase value		Plus key: increase value
Pic. 13: Decrease value		Minus key: decrease value
Pic. 14: Programming		Activate programming mode
Pic. 15: Program start/stop		Start/stop program
Pic. 16: Casting time		Start/stop casting time selector

#### 4.5.4 Display functions

	<ul style="list-style-type: none"> <li>In heating mode: Display of furnace temperature/holding time.</li> <li>In programming mode: Display of rate of increase/holding temperature/holding time.</li> </ul>
	Display of the actual program step (1 - 4)
	Display of the active program (1 - 9)
	<ul style="list-style-type: none"> <li>In Standby mode: Display of day (Mo = 1, Tue = 2, ...) and time (hh:mm).</li> <li>In heating mode: Display of time of completion</li> </ul>
	LED lights green for activated casting time selector.
	LED below the Start/Stop key lights green for activated heating operation.

#### 4.6 Delivery includes

The pre-heating furnace comes ready-to-use with:

- a ceramic insert,
- a PtRh-Pt thermo element and
- an exhaust pipe.

#### 4.7 Options

In place of the exhaust pipe the furnace may also be equipped with the following options for exhaust extraction and treatment:

- Mihm-Vogt vapour extractor DG2 or extraction hood DU1-DU3/3 with outdoor air outlet.
- Mihm-Vogt catalyst KN1 or KN2.

### 5. Transport, packing and storage

#### 5.1 Safety information for transport

##### Improper transport

##### CAUTION!

##### Damages due to improper transport!

In case of improper transport tangible damages in considerable extent may arise. Therefore:

- Use caution during offloading of the packages, during delivery and in-house transport and pay attention to the symbols and information on the packaging.
- Remove packaging only shortly before installation.

##### Weight of furnace

##### CAUTION! Risk of injuries due to improper lifting of the furnace!

Improper lifting of the furnace may cause back pain or injuries due to its high weight. Therefore:

- If necessary call additional person for lifting the furnace.

#### 5.2 Symbols on packaging

##### Top

The arrow points of the symbol mark the top of the package. They must always point to the top otherwise the contents could be damaged.

##### Fragile

Marks packages with breakable or sensitive contents. Handle the package with caution. Do not let drop and prevent from impacts.

#### 5.3 Transport inspection

Immediately after reception of the delivery check for completeness and transport damages.

In case of recognisable outside transport damages do the following:

- Do not accept delivery or only conditionally.
- Note scope of damages on the waybill or the delivery note of the forwarder.
- Initiate claim.

##### NOTE!

Claim any fault as soon as it becomes obvious. Claims for damages can only be submitted during the applicable claims period.

#### 5.4 Packaging

##### About packaging

The individual packages are packed according to the expected transport conditions. Only environmentally-friendly materials were used as packaging. Until assembly the packaging should protect the individual components from transport damages, corrosion and other damages. Therefore do not damage the packaging and remove only shortly before assembly.

##### Handling of packaging materials

Dispose of packaging materials according to the relevant legal regulations and local directives.


##### CAUTION! Environmental damages due to wrong disposal!

Packaging materials are valuable raw materials and may be further used in many cases or may be recycled and reused. Therefore:

- Dispose packaging in an environmentally-friendly way.
- Observe locally applicable disposal regulations. Possibly contract an authorised company for disposal.

## 5.5 Transport

### Transport of palettes with fork lifts

Pic. 21: Position of forks

Under the following conditions packages that are tightened on palettes may be transported with a forklift:

- The fork lift must be equipped according to the weight of the transport units.
- The driver must be authorised to drive the fork lift.

#### Loading:

1. Drive the forklift with the forks between or under the beams of the palette.
2. Drive the forks in until they show on the other side of the palette.
3. Make sure that the palette cannot tilt if the centre of gravity is off-centre.
4. Lift package and commence transport.

#### Storage of packages

Store packages under the following conditions:

- Do not keep in the open.
- Store in a dry and dust-free space.
- Do not expose to aggressive substances.
- Protect from solar radiation.
- Avoid mechanic vibration.

## 6. Installation and initial operation

### 6.1 Installation

#### 6.2 Safety

##### Personell

- Qualified laboratory worker, qualified dental technologist
- Certified electrician

#### General

##### **WARNING! Risk of injuries due to improper installation and initial operation!**

Improper installation and initial operation may lead to severe bodily injuries or tangible damages.

Therefore:

- Before work starts make sure that there is enough space for installation.
- Handle open sharp-edged parts with care.
- Take care of order and tidiness in the area of installation! Parts and tools that are lying around unattached are sources of accidents.
- Mount parts according to specifications.
- Secure parts to prevent them from falling or toppling over.

#### Electric system

##### **DANGER! Risk of fatal injury from electrical current!**

The imminent risk of fatal injury exists if live parts are touched. Therefore:

- Before any work switch off mains supply and secure against reconnection.

#### Weight of furnace

##### **CAUTION! Risk of injuries due to improper lifting of the furnace!**

Due to its high weight improper lifting of the furnace may cause back pain and injuries. Therefore:

- Exclusively lift the furnace as described in section "Transport, packaging and storage".
- If necessary call an additional person for lifting the furnace.

## 6.3 Installation

### 6.3.1 Installation of pre-heating furnace

#### Installation of furnace

Install the pre-heating furnace on an even surface. Make sure that the installation surface of the pre-heating furnace has sufficient carrying capacity for the weight of the pre-heating furnace. Observe the conditions of operation (» section Technical data).

##### **CAUTION! Risk of ignition for objects in the vicinity due to high unit temperature!**

Even during proper use parts of the pre-heating furnace may have high temperatures. Combustible, highly flammable or explosive materials in the vicinity might be ignited. Therefore:

- Make sure that furniture and fittings are not made of highly flammable or explosive materials.
- Do not store any highly flammable or combustible gases or liquids in the installation area of the pre-heating furnace

#### Extraction of vapours/gases

To exclude an accumulation of vapours and gases in the room of installation equip the laboratory furnace with a vapour extractor or catalyst and install under an extraction hood (» section Installation of vapour extractor and catalyst, connection of extraction hood).

## 6.4 Connection to mains supply

#### Power connection

- 230 V-earthed socket

Thereby check for a secure protective earth conductor connection! The earth conductor contacts in the earthed socket must not be contaminated. Observe requirements for electric connections (» connection values in section Technical data).

## 6.5.4 Adjustment of furnace parameter

The pre-heating furnace is delivered with

- pre-set time
- with standard parameters.

At initial operation it might be necessary to adjust time and other furnace parameters to the individual requirements.

**NOTE!** The FINOFORNO controller does not automatically adjust to daylight saving time.

#### Select and quit parameter mode

1. Keep **Prg** Key pressed.
2. Switch on furnaces using main switch.

3. Release **Prg** key.

Parameter mode is active.

Select parameters with **↓** and **↗** keys.

Modify parameters with **Plus** and **Minus** keys.

#### Parameter Adjustment

P1	Time mode: 24h-mode/12h-mode
P2	Time: day (1: Monday, 2: Tuesday, ...)
P4	Time: hour
P5	Time: minute
P6	Beep: on (1, standard)/off (0)
P7	Temperature unit: °C/°F

In case of connected vapour extractor, catalyst or extraction hood:

- |    |   |
|----|---|
| P8 | Switch-off temperature for socket "Gebläse" |
| P9 | Switch-off temperature for socket "Dunst."  |

**NOTE!**

If Speed investment materials are used it is recommended to select the switch-off temperature approx. 20 °C higher than the final temperature indicated by the manufacturer. This ensures continuous operation of the units for furnace exhaust extraction.

#### Quitting parameter mode

1. Select parameter P8.
2. Press **↗** key.

On quitting the parameter mode parameter modifications are automatically saved.

## 6.5.5 Run revision program

For the calculation of the program start time (» section Programming) the controller assumes an average voltage of 230 V. Due to local circumstances the actual voltage value may be different. This may lead to faulty time calculations of the controller. The revision program compensates these deviations. If you notice larger time differences during operation of the pre-heating furnace (> 30 minutes) run revision program as follows:

**NOTE!** The running time of the revision program is approx. 8 hours.

1. Keep Start/Stop key pressed.
2. Switch on furnaces using main switch.
3. Release Start/Stop key.

## 7. Operation

### 7.1 Safety

#### Personnel

- Qualified laboratory worker,  
qualified dental technologist

#### Basics

**WARNING!** Risk of injuries due to improper operation!

Improper operation may lead to severe bodily injuries or tangible damages. Therefore:

- Follow all operational steps according to the specifications in this user manual.
- Before any work make sure that all covers and safety devices are installed and in orderly working condition.
- Do never put any safety devices out of service during operation.

**CAUTION!** Risk of burns due to hot surfaces at the furnace door!

During normal operation the surface temperature at the furnace door may exceed 60°C. Even a switched off furnace may have a high residual heat from a previous heating process.

Therefore:

- Use sufficiently heat-insulated protective gloves.
- Open the furnace door only with therefore designed lever/handle.

#### Heating chamber

**CAUTION! Risk of burns due to hot surfaces in the heating chamber!**

Even a switched off furnace may have a high residual heat from a previous heating process. Therefore:

- To load or remove fired objects strictly use tongs that are sufficiently long.
- Use sufficiently heat-insulated protective gloves.

**CAUTION! Risk of damage for door inserts due to improper handling!**

The ceramic of the door insert is highly porous and therefore sensitive to scratches and impacts.

Therefore:

- Always treat the door insert carefully.
- Avoid touching with zircon block and tongs.

**CAUTION! Risk of damage of the controller electronics due to electrical discharge!**

During heating operation there is the risk that if the distance between metallic objects and the walls of the heating chamber is too small an electrical discharge may take place. This may damage the controller electronics of the furnace beyond repair. Therefore:

- Make sure that metallic objects do not touch the walls of the heating chamber.

## 7.2 Loading of furnace

	<ol style="list-style-type: none"> <li>1. Open furnace door.</li> <li>2. Place ceramic insert on floor of heating chamber.</li> </ol> <p><b>NOTE!</b> In any case use a ceramic insert to protect the floor of the heating chamber against contamination.</p> <ol style="list-style-type: none"> <li>3. Distribute casting moulds in the heating chamber in a way that they do not touch the walls of the heating chamber.</li> <li>4. Close furnace door.</li> </ol>
---	---

Abb. 24: Ceramic insert

## 7.3 Switch on furnace

	<ul style="list-style-type: none"> <li>• Press main switch. The mains control light in the main switch is on.</li> </ul> <p>After a few seconds the operating display shows the furnace temperature (here e.g.: 27 °C).</p>
	<p>Pic. 25: Furnace temperature</p> <p>The day/time display is alternating between time and weekday (her e.g.: Wednesday, 08:32 hours).</p>
	<p>Pic. 26: Display time</p> <p>Pic. 27: Display weekday</p>

## 7.4 Programming

<b>Activate programming mode, select program</b>	The controller offers the option to store up to nine different heating programs.
 Pic. 28: LED Prog.	1. Press $\blacktriangleright$ key until the LED below "Program number" display is on. 2. Select program (1-9) with <b>Plus</b> and <b>Minus</b> keys (here e.g.: program 1).
 Pic. 29: Operating display, LED blinks	3. Keep the <b>Prg</b> key pressed until the operating display indicates ready-for-input. The LED blinks. <b>NOTE!</b> If no selection has been made within 10 seconds using the Prg key the displays switches back to the normal display mode.
<b>Programming of steps, temperatures, holding times</b>	In the programming mode individual heating curves can be compiled by the selection of up to 4 steps.
 Pic. 30: Step 1	Enter heating parameter for step 1:
 Pic. 31: Enter heating speed	1. Select heating speed (C/min) with <b>Plus</b> and <b>Minus</b> keys (here e.g.: 9 °C/min). <b>NOTE!</b> If "0" is selected the furnace heats with the maximum heating speed (> 9 °C/min) non-linear.
 Pic. 32: Enter temperature	2. Select temperature (C) with $\blacktriangleright$ . 3. Enter temperature with <b>Plus</b> and <b>Minus</b> keys (here e.g.: 1150 C). <b>NOTE!</b> Maximum temperature: 1500 °C.
 Pic. 33: Enter holding time	4. Select holding time (min) with $\blacktriangleright$ key. 5. Enter holding time with <b>Plus</b> and <b>Minus</b> keys. <b>NOTE!</b> Maximum holding time: 99 minutes. After holding time has terminated you will hear a beep. Nevertheless the furnace will hold the temperature until the operator quits the program.
 Pic. 34: Additional steps	6. If necessary select another step using $\blacktriangleright$ key and <b>Plus</b> and <b>Minus</b> keys (here e.g. second step). 7. Repeat program steps 1 to 5 for more steps.
<b>Quit programming mode</b>	After the respective heating parameters for a step have been programmed the programming mode can be quitted. To quit the programming mode keep the <b>Prg</b> key pressed until the LED below the operation display is permanently on. The program is now permanently saved. <b>NOTE!</b> If a program should only consist of one single step it is important to make sure that all values in the further steps are set to "0".

## 7.5 Select program

	<ol style="list-style-type: none"> <li>Press <math>\triangleright</math> key until the LED below the "Prog" display is on.</li> <li>Select program (1-9) with the <b>Plus</b> and <b>Minus</b> keys.</li> </ol>
Pic. 35: LED Prog.	

## 7.6 Start program

	<ul style="list-style-type: none"> <li>Start program with <b>Start/Stop</b> key. The program starts and the furnace start heating.</li> <li>During the program cycle the operation display shows the furnace temperature respectively the remaining holding time.</li> <li>The day/time display shows the time of completion of the program.</li> <li>Using the <math>\triangleright</math> key you can display the set temperature/holding time.</li> </ul> <p>At any time the program can be aborted with the <b>Start/Stop</b> key. After the <b>Start/Stop</b> key has been pressed again the program will start from the beginning.</p>
--	--

## 7.7 Start program automatically

	<p>The furnace can be started automatically with an integrated timer in a way that the heating program is finished at a pre-selected completion time. With the timer you can determine the completion time.</p>
	<ol style="list-style-type: none"> <li>Select program (» subsection select program).</li> </ol>

Pic. 36: Select program

	<ol style="list-style-type: none"> <li>Activate timer with <math>\odot</math> key. The green LED below the <math>\odot</math> key lights up. The yellow LED "week day" blinks.</li> </ol>
Pic. 34: Timer activated	

	<ol style="list-style-type: none"> <li>Select week day (Mo = 1, Tue = 2, ...) with <b>Plus</b> and <b>Minus</b> keys (here example Wednesday).</li> </ol>
Pic. 35: Select week day	

	<ol style="list-style-type: none"> <li>Press <math>\triangleright</math> Key. The yellow LED "Time" blinks.</li> <li>Select time (hours) with <b>Plus</b> and <b>Minus</b> keys (here example six o'clock).</li> </ol>
Pic. 36: Select hour	

	<ol style="list-style-type: none"> <li>Press <math>\triangleright</math> Key.</li> <li>Select time (minutes) with <b>Plus</b> and <b>Minus</b> keys (here example 30 minutes).</li> </ol>
Pic. 37: Select minute	

	<ol style="list-style-type: none"> <li>Press <math>\triangleright</math> key to finish entry of the completion time.</li> </ol>
Pic. 38: Casting time	

	<p>After the casting time has been entered the start time can be checked by pressing the <b>Start/Stop</b> key.</p> <p> <b>NOTE!</b> By pressing the <math>\odot</math> key again the integrated timer is deactivated. This makes the immediate manual start of the selected program possible.</p>
Pic. 39: Check start time	

## 7.8 Temperatur stepping

Casting moulds with different final temperatures can be heated in the furnace successively (final temperature A and higher final temperature B).

If a program should only consist of one step (e.g. for „Shock“ heating) make absolutely sure that all values in further steps are set to "0"!

1. Compile heating program where final temperature A is reached in the last step (» section Programming).
2. Remove casting moulds with final temperature A from the heating chamber after the heating program has terminated.
3. Keep Prg key pressed.  
Now the programming mode is activated.
4. Increase the temperature of the last step to temperature B with the Plus key.
5. Keep Prg key pressed to quit the programming mode. The heating program repeats the last step with the new temperature B.

6. After heating program is finished remove casting moulds with final temperature B from the heating chamber.

## 7.9 Switch off

- To switch off the furnace press the main switch. The mains control light in the main switch is off.

**CAUTION! Risk of burns due to residual heat of the pre-heating furnace!**

Even if the furnace is switched off the heating chamber may still have a considerable residual heat. There is the risk of burns at the walls of the heating chamber and at the furnace door. Therefore:

- Before working on the furnace make sure that it has cooled sufficiently. The furnace needs at least 12 hours to cool from maximum temperature to room temperature.

## 8. Malfunctions

In the following section possible causes for malfunctions and processes for their repair are described. In case that malfunctions cannot be eliminated referring to the following information please contact the manufacturer (» Service address, last page).

## 8.1 Safety

### Personell

- Some of the here described processes for troubleshooting may, if not marked differently, be performed by the operator.
- Work on the electrical system may exclusively be performed by certified electricians.
- Some repairs may only be done by the service department of FINO. In the description of individual malfunctions these are separately indicated.

### Personal protective equipment

#### NOTE!

In the warning notes of this section protective equipment that has to be worn for specific work is specifically indicated.

### Electrical system

#### DANGER! Risk of fatal injury from electrical current!

On contact with live parts there is the risk fatal injury.

Therefore:

- Before any work switch off the mains supply and pull the mains plug.

#### CAUTION! Risk of burns due to residual heat of the pre-heating furnace!

Even if the furnace is switched off the heating chamber may still have a considerable residual heat. There is the risk of burns at the walls of the heating chamber and at the furnace door. Therefore:

- Before working on the furnace make sure that it has cooled sufficiently. The furnace needs at least 12 hours to cool from maximum temperature to room temperature.

## 8.2 Malfunction indications

The following devices indicate malfunctions:

### Operating display

The operating display shows error messages in case the electronics are malfunctioning. These are described in the following table of malfunctions.

## 8.3 Table of malfunctions

Malfunction	Possible cause	Troubleshooting	Performed by
Wrong time	Time is wrongly stored by the controller	Set time correctly (» section Initial operation).	Operator
No display, power light on	Defunct fuse	Switch furnace off, wait for 30 seconds, switch on again. In case of repeated malfunction let controller be exchanged.	
No display, power light off	No mains power available	Check on site circuit breakers and power connection. If necessary call certified electrician.	
Furnace does not start	Powercut/interruption of power supply	Check if mains connection is interrupted.	

Malfunction	Possible cause	Troubleshooting	Performed by
Display "Er01"	Defect thermo element	Exchange thermo element.	Certified electrician, Service
Display "Er02"	Loose thermo element connections	Retighten thermo element connections.	
Display "Er04"	Wrong connection/reverse polarity of thermo element	Change thermo element connections.	
Display "Er05"	Thermo element shortcut or heating is defect	Check heating/thermo element and exchange if necessary.	
Display "Er03"	Defunct amplifier electronics	Controller must be exchanged.	Service
Display "Er06"	Electronics defect	Check electronics and exchange if necessary.	
Display "--:--"	Door safety switch is jammed or defect	Visual check of door switch. Call Service.	
Display "Furnace is heating", but furnace does not	Defunct heating	Conduction of heating must be checked.	
	Defunct controller	Controller must be exchanged.	
Heating programs and programs are not permanently saved	Defunct controller	Controller must be exchanged.	

## **8.4 System information for malfunction diagnosis Carried out by the operator**

In case of malfunctions of the controller electronics and irregularities in the controller software it is necessary to contact the service (» Table of malfunctions).

For a closer malfunction diagnosis the service department needs precise system information will ask them from the user:

- Furnace type
- Serial number
- Hardware version
- Software version

Retrieve system information as follows:

1. Switch furnace on with mains switch.  
For approx. 3 seconds numbers and/or text will appear in the display.
2. Write this information down and pass on to the service department.

## **9. Maintenance**

### **9.1 Safety**

#### **Personell**

- The here described maintenance can be done by the user.

#### **Residual heat**

##### **CAUTION! Risk of burns due to residual heat of the pre-heating furnace!**

Even if the furnace is switched off the heating chamber may still have a considerable residual heat. There is the risk of burns at the walls of the heating chamber and at the furnace door. Therefore:

- Before working on the furnace make sure that it has cooled sufficiently. The furnace needs at least 12 hours to cool from maximum temperature to room temperature.

#### **Electrical system**

##### **DANGER! Risk of fatal injury from electrical current!**

On contact with live parts there is the risk fatal injury.

Therefore:

- Before any work switch off the mains supply and pull the mains plug.

## **9.2 Maintenance schedule**

In the following sections maintenance work is described that is required for an optimum and problem-free operation. In case that during regular checks an increased wear is observed shorten the required maintenance intervals according to the respective actual wear and tear. If you should have questions regarding maintenance work and intervals contact the manufacturer, see Service address, last page.

**Maintenance:** clean heating chamber

**Interval:** daily

**To be performed by:** Operator

**Maintenance:** heat empty furnace

**Interval:** in the first three month of operation every fortnight

**To be performed by:** Operator

**Maintenance:** heat empty furnace

**Interval:** After the initial month of operation every month

**To be performed by:** Operator

## **9.3 Maintenance work**

### **9.3.1 Cleaning of furnace**

Keep the heating chamber always clean. Regularly remove slag residues and other foreign bodies from the heating chamber with a soft brush.

### **9.3.2 Heating of empty furnace**

Heat the unloaded furnace for 90 minutes to 1050 °C.

**NOTE!**  
This process maintains the protective oxide layer on the heating wire. Additionally potential wax residues are removed from the heating chamber.

## **10. Disassembling and disposal**

After the end of the service-life disassemble the unit and dispose in an environmentally-orientated way.

### **10.1 Safety**

#### **Personell**

- The here described maintenance can be done by the user.

## **10.4 Disposal according to EC directive 2002/96/EC**

 This unit must not be disposed of with household waste! Please contact the manufacturer or your supplier!

## **10.5 Registration number**

WEEE-Reg.-No. DE 75542610

#### **Residual heat**

##### **CAUTION! Risk of burns due to residual heat of the pre-heating furnace!**

Even if the furnace is switched off the heating chamber may still have a considerable residual heat. There is the risk of burns at the walls of the heating chamber and at the furnace door. Therefore:

- Before working on the furnace make sure that it has cooled sufficiently. The furnace needs at least 12 hours to cool from maximum temperature to room temperature.

## **10.2 Disassembling**

About the disassembling of the pre-heating furnace:

- Switch pre-heating furnace off and disconnect from mains.
- Remove all connections as far as they are installed from the pre-heating furnace:
  - Vapour extractor/catalyst
  - Extractor hood
  - Extractor pipe

## **10.3 Disposal**

The pre-heating furnace including options is classified as a medical product according to the electric devices directive and must not be disposed of with household waste.

##### **CAUTION! Risk of environmental pollution due to wrong disposal!**

All components of the pre-heating furnace are subject to special refuse treatment and may only be disposed by authorised companies.

The local municipal authorities or specific certified disposal companies will inform about environmentally friendly disposal according to regulations.

FINO GmbH  
Mangelsfeld 18  
D-97708 Bad Bocklet  
Tel +49-97 08-90 94 20  
Fax +49-97 08-90 94 21  
info@fino.com · www.fino.com

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