FINO

IMPULS – CONNECTOR

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Deutsch

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9.1 Registrierungsnummer
WEEE-Reg.-Nr. DE 75542610

FINO
IMPULS – CONNECTOR
Spot-Welding Unit 72490

English

User Information

1. Fields of application

The FINO IMPULS - CONNECTOR is used for spot welding respectively for the making of welding joints of all dental alloys.

2. Advantages

- Suitable for all dental alloys
- High welding performance
- Short charging times
- Absolutely maintenance-free

3. Safety precautions

It is recommended to wear the included goggles whenever you work with the FINO IMPULS - CONNECTOR.

There is the risk of burns if you touch fresh welding joints.

4. Processing instructions

This unit is very powerful (88.000 MF). The achieved high welding performance guarantees stable welding joints for all dental alloys. Additionally the fine adjustment enables welding for different applications from a broken bar to orthodontic wires. Fast welding sequences are possible due to the extremely short charging times. Time-consuming fixations with resin or soldering investment material are not longer required.

4.1 Getting started

After unpacking of the unit it should immediately be checked for transport damages.

If there are any transport damages these should be claimed immediately. The unit is delivered ready-to-use. After making sure that the voltage on the type label of the unit corresponds to the available mains supply the unit plug is connected to the socket.

4.2 Operation

The switch is put on position "ON" so that the green control light in the mains switch lights up.

The three flexible cables are plugged into the unit's jacks. The two black ones in the black jacks and the red cable into the red jack. To fixate a part directly to a workpiece the black cable is connected with the holding clamp to the base (workpiece). With the clip of the red cable the part that must be fastened is taken up and both parts are merged to each other. This process is done in switch position „Single arrow“.

To fuse 3 parts with each other or 2 parts with a large connecting area where a third part e.g. a wire pin that is used as a connecting joint the two black cables are connected to one base respectively and the red cable holds the connecting part. This process is done in switch position „Double arrow“. In this case the control electronics splits the energy output in two equal parts for each black cable.

Thereby it is ensured that a spot welding joint irrespective of connecting area the connecting pressure is equally created on the left as well as the right part. The rocker switch with the single and double arrow should be positioned depending on the alloy that must be joined.

Metal and alloys like silver and gold have the highest electrical conductivity and require therefore the respective low specific resistance. A possible connection "gold to gold" or "silver to silver" should be done in switch position „Double arrow“.

Thereby the energy is released in two steps and the fusion of the metal surfaces is more intensive. Due to the high conductivity the spot-welding effect is reduced, which means an electric arc only forms to a lower degree. In this case a double impulse counteracts this. Other metals with reduced conductivity or connections "gold to steel" can be done in switch position „Single arrow“.

The impulse strength can be variably adjusted with the regulator e.g. fixing of a wire pin with:

Ø 0,5 mm-Power approx. 15-20 Volt
Ø 0,8 mm-Power approx. 20-25 Volt
Ø 1,2 mm-Power approx. 25-35 Volt

The larger the connecting area of the two parts is the higher the impulse strength must be selected. The ideal distance of an applied holding clip to the welding point is approx. 10-20 mm. If the parts that are to be fixed are held together the pushing of the foot switch releases an electric arc. After the discharge of the energy storage in the unit the red light indicates that the unit is automatically recharging. If the regulator is e.g. set to the maximum value of 40 Volt the charging time is only approx. 10 seconds until the green light again indicates „Ready“. For delicate, precise small parts the medium respectively small adapter clips can be used.

4.3 Maintenance

The unit is maintenance-free. Repairs of the unit should only be done by FINO or an authorised workshop. The user is liable for damages that are caused by arbitrary repairs. The terms of delivery and guarantee of DT&Shop are unrestricted applicable for this unit.

5. Storage

No special measures required.

6. Technical data

Voltage: 220 – 240V, 50 Hz
Secondary voltage: max. 40V
Storage capacity: 88.000 MF
Exterior dimensions: 85 x 180 x 235 mm
Weight: 3,4 kg

7. Delivery forms

FINO IMPULS- CONNECTOR

Spot-Welding Unit each 72490

7.1 Accessories

FINO Impulse C. Basic cable set
each 72491
FINO Impulse C. Crocodile clip (large)
each 72492
FINO Impulse C. Copper tip
each 72493
Adapter cable with medium clips
each 72494
Adapter cable with small clips
each 72495

8. Guarantee

Our technical recommendations of application are based on our own experiences and tests and should only be regarded as guidelines. It rests with the skills and experience of the user to verify that the products supplied by us are suitable for the intended procedures. Our products are undergoing a continuous further development. We reserve the right of changes in construction and composition. It is understood that we guarantee the impeccable quality of our products.

9. Disposal in accordance with EU guideline 2002/96/EU

This equipment may not be disposed in domestic waste! Please contact the manufacturer or your dealer!

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