Stereomikroskope

MSL4000-Serie

KSW4000-Serie

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Stereo Microscopes

MSL4000 Series

KSW4000 Series

Operating Manual

Inhalt

MSL4000-Serie

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Index

MSL4000 Series

Introduction and equipment summary	25
1. Description	26
2. Unpacking and assembly	27
3. Supporting plate	27
4. Focusing	27
5. Adjustment of eye distance and diopter	28
6. Illumination	28
7. Change of magnification	29
8. Technical data	29

KSW4000 Series

Introduction and equipment summary	31
1. Description	32
2. Assembly vertical	33
3. Operation	34
4. Assembly horizontal	35
5. Bulb replacement	37
6. Technical data	38
Maintenance	39
Reclamation and recycling	39
Warranty	40
Warranty extension	41

The MSL4000 Series

Congratulations on your purchase of a new KRÜSS microscope! You now own a device of the highest optical as well as mechanical quality. The excellent operating convenience of our microscopes will guarantee many years of reliable and trouble-free operation. Please follow the instructions on page 19 in chapter *Maintenance*.


A. KRÜSS would like to point out that this instruction manual contains important safety and maintenance information and should therefore be made available to all users. A.KRÜSS does not accept any liability for the improper use of the microscope.

The MSL4000 Series includes the following microscopes:

MSL4000-10/30-IL-TL	Magnification: 10x and 30x Illumination: LED incident light and transmitted light
MSL4000-10/30-IL-S	Magnification: 10x and 30x Illumination: Light bulb 12 V 20 W Swivel arm
MSL4000-10/30-S	Magnification: 10x and 30x Swivel arm
MSL4000-20/40-IL-TL	Magnification: 20x and 40x Illumination: LED incident light and transmitted light
MSL4000-20/40-IL-S	Magnification: 20x and 40x Illumination: Light bulb 12 V 20 W Swivel arm
MSL4000-20/40-S	Magnification: 20x and 40x Swivel arm

1. Description MSL4000

2. Unpacking and assembly

The Styrofoam case contains the following parts:

- Microscope stand with microscope head
- Eye pieces and lenses attached to the microscope head
- B/W contrast plate as stand base (only for models without swivel arm)
- Transparent glass plate (only for models with transmitted light and without swivelling arm)
- Vinyl protecting cap
- Instruction manual

Remove the microscope from the Styrofoam case. The eyepieces are already in place.

3. Supporting plate

The supporting plate is to be selected according to the material to be examined. The transparent glass plate is used for transparent materials, the white plate for coloured and the black plate for white materials.

4. Focusing

- Place the material to be examined on the supporting plate and illuminate, if necessary.
- While looking through the microscope, release the locking screw and move the head up and down until the image is in focus. Now, fasten the screw.

You can get a sharp image of the material by slowly turning the focus adjustment knob.

5. Adjustment of eye distance and diopter

- Move both eyepieces until they match your eye distance and until you can see an observation field through the microscope.

The observer might have small diopter differences of the eyes which have to be compensated in order to avoid eye fatigue.

- Look with your right eye through the right eyepiece and focus the image using the focus adjustment knob.
- Now look with your left eye through the left eyepiece and focus the image using the diopter adjustment without turning the focus adjustment knob.

6. Illumination

Some microscope models have been fitted with an illumination device that is simply turned on by a switch located on the housing.

- The main switch switches all of the illumination on or off.
- The 3-way switch in the neutral position is the OFF position. By moving the switch forward, you turn on the incident light, by moving backwards, you turn on the transmitted light.

7. Change of magnification

The microscopes of the MSL4000 series have been fitted with a turret lens to select two magnifications. The magnification is changed by rotating the lens housing by 180°. The respective magnification is then displayed. The image should first be focused with a low magnification and you should then switch to a higher magnification.

8. Technical data

General

- Adjustable eye distance
- 45° inclined eyepiece
- Metal stand

Optical features

- Eyepieces: 10x wide field
20x wide field (optional)
- Object lenses: 1x and 3x or
2x and 4x

Illumination

- 6 V LED incident and transmitted light or incident light with swivel arm
12 V 20 W light bulb

Power supply

- Universal power transformer 90 – 240 V, 50 Hz
Output: 6 V DC

The KSW4000 Series

Congratulations on your purchase of a new KRÜSS microscope! You now own a device of the highest optical as well as mechanical quality. The excellent operating convenience of our microscopes will guarantee many years of reliable and trouble-free operation. Please follow the instructions on page 19 in chapter *Maintenance*.

A. KRÜSS would like to point out that this instruction manual contains important safety and maintenance information and should therefore be made available to all users. A.KRÜSS does not accept any liability for the improper use of the microscope.

The KSW4000 Series includes the following microscopes:

KSW4000

Magnification: 10x and 30x

Illumination: LED incident light and transmitted light

KSW4000-K

Magnification: 10x and 30x

Illumination: 12 V 20 W light bulb (transmitted light)
12 V 20 W halogen bulb (incident light)


KSW4000-K-W

Magnification: 10x and 30x

Illumination: 12 V 20 W light bulb (transmitted light)
12 V 20 W halogen bulb (incident light)

Horizontal assembly is possible

1. Description KSW4000

2. Assembly vertical

- Attach the cups to the eyepieces.

- Insert the white plate for coloured, the black plate for white materials.

- Insert the dark field condenser for the examination of precious stones in the dark field.

- Insert stone tong into socket.

The KSW4000 has been fitted with 2 sockets in order to be used by both right-handed and left-handed people.

- For incident light, insert the lightguide into the lamp holder on the rear side.

- Insert the plug into the mains connection jack and the transformer into the 220 V outlet.

3. Operation

- Place the material to be examined on the supporting plate and illuminate, if applicable.
- Look through the eyepiece.
- Release the locking screw of the head and move the head up and down until there is a clear image. Now fasten the screw again.
- Turn the focus adjustment knob slowly, to get a sharp image of the material.
- Move both eyepieces until they match your eye distance.
- Look with your right eye through the right eyepiece and focus the image using the focus adjustment knob.
- Now look with your left eye through the left eyepiece and focus the image using the diopter adjustment without changing the focus adjustment.
- The main switch turns all of the illumination on or off. The 3-way switch in the neutral position is the OFF position. By moving the switch forward, you turn on the incident light, by moving backwards, you turn on the transmitted light.

- The magnification is changed by rotating the lens housing by 180°. The respective magnification is then displayed. The image should first be focused with a low magnification and you should then switch to a higher magnification.

4. Assembly horizontal

- Insert the pillar into the block horizontally. Fasten by using the screw on the side.

- Insert the cell stage vertically into the rear hole of the block. Fasten by using the screw.

- Attach the polarization device to the cell stage.

- Insert the stone tong into the socket of the cell stage.

- Place optical cell on the stage.

- Insert the plug into the mains connection jack and the transformer into the 220 V outlet.

5. Bulb replacement

Important note!

Always pull the mains plug before replacing the bulb.

Important note!

When touching the bulb, use a soft cloth.
Finger prints on the glass bulb affect the light output and reduce the service life.

Transmitted light

- Remove the supporting plate
- Remove the defective bulb from the socket
- Insert a new bulb into the socket
- Insert the supporting plate

Incident light (cold light source)

- Loosen the two screws
- Remove the lamp holder

- Pull off the socket
- Release the bayonet catch

- Remove the defective bulb
- Insert the new bulb between iris diaphragm and spring
- Fasten by using the bayonet catch
- Attach the bulb socket
- Insert the lamp holder and fasten by means of the two screws

6. Technical data

General

- Adjustable eye distance
- 45° inclined eyepiece
- Metal stand

Optical features

- Eyepieces: 10x wide field
20x wide field (optionally)
- Object lenses: 1x and 3x

Illumination

- KSW4000: 6 V LED incident and transmitted light
KSW4000-K and -K-W: 12 V 20 W transmitted light
12 V 20 W cold light source

Power supply

- Universal mains transformer 90–240 V, 50 Hz
KSW4000 output: 6 V DC
KSW4000-K and -K-W output: 12 V DC

Maintenance

Just like all precision instruments, your microscope should be handled with care. Protect your microscope from dust, condensation, vibration and direct exposure to sunlight.

Keep the optical system very clean. Remove dust with a soft brush or clean air.

Remove finger prints with a lint-free cloth soaked in an ether/alcohol mixture (80% ether and 20% alcohol). Special cleaning paper can be purchased in a photo store.

Do not use any corrosive solutions for the microscope!

After usage, place the microscope back into the box and protect it with a protective dust cover.

Important note!

Aggressive substances may damage the finish or the plastic casing!

Important note!

In case of a malfunction, the microscope should be repaired by A. KRÜSS Optronic or an authorized repair shop

Reclamation and recycling

The device could be an important source of raw material. Please do not dispose as garbage but collect it separately for recycling and reclamation of the contained materials. In case of improper disposal, the materials might be harmful for the environment and health.

The manufacturer of the microscope, A. KRÜSS Optronic GmbH collects, uses and recycles the contained raw materials. However, the reclamation requires your support.

Should you decide to dispose this microscope, do not attempt to open it or to use parts of it for other purposes than described in this instruction manual. Take the device to the dealer from which you have bought it. The dealer should accept the microscope free of charge.

The reclamation of the raw materials takes place in compliance with the European Directive 2002/96/EC and all other applicable requirements.

Warranty

A. KRÜSS Optronic GmbH accepts the guarantee for material and manufacture of the microscope for a period of 24 months from the date of shipment. During this term of warranty A. KRÜSS Optronic GmbH shall remedy any defect by repairing or replacing the device if subject to the warranty claim.

For repair or service, the device must be sent to A. KRÜSS Optronic. In case of warranty repairs, A. KRÜSS Optronic GmbH bears the costs of the shipment; otherwise it is at the expense of the customer.

A. KRÜSS Optronic guarantees that the hardware designed by A. KRÜSS for this device is free from defects when used according to the information provided by the manufacturer.

However, A. KRÜSS Optronic does not guarantee the correct and uninterrupted operation of the device or that this instruction manual is error-free.

We are also not liable for subsequent damage.

Warranty limitations:

The warranty described above does not include errors and defects resulting from improper use, modifications, misuse or any operation outside the specified environment or unauthorized maintenance.

Further claims will not be accepted. A. KRÜSS Optronic will explicitly not guarantee the application or efficient use for certain applications.

A. KRÜSS reserves the right to modify this instruction manual and the technical data of the described device at any time.

KRÜSS microscopes are only fit for shipment if properly packed into the complete original packaging. If needed, request a replacement packaging from your supplier.

A.KRÜSS Optronic GmbH
Alsterdorfer Strasse 276–278
22297 Hamburg / Germany
Tel. +49-(0) 40-514317-0
Fax. +49-(0) 40-514317-60
Internet: www.kruess.com
Email: service@kruess.com

We value your opinion!

We grant you a 1 year warranty extension period if you return the completely filled out questionnaire to us!

Company _____
 Contact person _____
 Email _____
 Address _____
 Serial number _____

Was the date of delivery kept? yes no

Did any problems occur during start-up? yes no

Are you already a customer a new customer

Are you a distributor an end customer a university/institute

How did you get to know about us? recommendation _____
 catalogue _____
 magazine _____
 exhibition _____
 internet/other _____

Which laboratory magazines do you read? _____

Which catalogues do you read? _____

Would you like to be informed about our latest updates? yes no

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A.KRÜSS Optronik GmbH