## Materialbeschreibung

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Rx only

## Description of material

Laborsil pro 70 / 90 are 1:1 addition-vulcanising kneadable silicones for the universal use in the dental laboratory with a hardness of 70 and 90 Shore A.

## Indication

For the production of copings, of repair models with under-growing areas, to block out undercuts and embedding silicone for the denture pressure procedure.

## Technical data according to DIN EN ISO 4823

Consistency	kneadable – Typ 0
Working time	1:15 min.
Demoulding at 23 °C	> 6:00 min.
Deformation under pressure	3,0%
Elastic recovery	> 99,2 %
Linear dimensional change	< 0,2 %
Final hardness	70 / 90 Shore A

The mixing and processing times refer to a room temperature of 23 °C and a relative air humidity of 50 %. Lower temperatures prolong, higher temperatures shorten these times.

## Processing mixing

Take two equal-sized amounts of components A and B and knead them manually until a homogenous colour develops.

## Important advice

Polymerisation may be disturbed by use of e.g. latex gloves, creams, cleaning agents, resins etc.

## Keep away from children!

**Only for use in the dental laboratory!**

## Storage condition

At room temperature, i. e. 18 °C–28 °C

**Durability** see date of expiry.

## Forms of delivery

### Laborsil pro 70

 REF	<b>D38390</b>	2 x 750 g	tub
 REF	<b>D38391</b>	2 x 1.9 kg	tub
 REF	<b>D38392</b>	2 x 5 kg	bucket

### Laborsil pro 90

 REF	<b>D38393</b>	2 x 750 g	tub
 REF	<b>D38394</b>	2 x 1.9 kg	tub
 REF	<b>D38395</b>	2 x 5 kg	bucket

This product was developed for dental lab use and has to be applied according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. The user is generally obligated to check the product before use. Especially if the purpose is not mentioned in this instructions.

Date of information: 2020.02.13