

Server instruction

IP polish motor PM 30N

General instructions

Before starting, maintenance- or service works of the machine, it is advisable, to read the present handout very carefully.

The safety of the user as well as the present persons depends on a conscientious and carefully behaviour.

By using the machine it is very important, to know everything about the position and function of all steering preparation.

- the server must have a good experience about the above mentioned machine. All parts of the machine ,which are not strong in wear and tear, should be checked periodical.
- never damage the electrical installation.
- If there is some working interruption ,please make sure, that all steering preparations are in the neutral position.
- The tool has a safety system, which should never be changed.
- The enable and maintenance works should only made by authorized persons.
- For cleaning from the machine parts do not use petrol, solution means or any light inflammable means: please only use no poisonous and no inflammable means.

Rest danger

If you will follow to the general safety rules ,the machines of the company Dental IP will have no other, rest dangers, either these, which are caused by the user himself.

Dispatch of the correspondence

If you want to ask the producer for some information,please specify the following datas:

- machine model
- serial number
- year of building
- purchase date
- exactly information about the find out mistake
- time of working,numbers of the working hours

Basis

The aim of the present handout is, to ease the relationship between the client and the after – sale – service of the company IP Division GmbH and particularly to publish the server instruction of the machine.

Every machine will be delivered with the belonging to it handout ,the in it included service instruction represents the necessary gurantee condition for a correct working of the machine.

It is therefore very important to support the following rules very conscientious.

Attention: I fy you do not take notice of the rules in the present handout, you will have no more any guarantee.

Transport and unpacking

Please unpack the parts in the box very carrefully You should check all parts, for discover possible damages ,which can happen by the dispatch.

If there are some damages you should inform immediately the producer about it.

It is very important to specify the exact name of the damaged piece.

Working definitions

- connect the powercable onto the wall socket and control, whether the supply tension will be the same.
- Attention: do not work with free material and/or free article of clothings.
- If you do not observe the rules, there will be a rest risk, which will twist a little bit.
- In screw the delivered polish spindle . In screw right or left hand side.
- Start the polishmachine with the help of the switch on the foot.

Technical datas:

Velocity: 1400/2800/U/m

Weight: 10.00 Kg

Maintenance: 220 V 50/60 Hz

Power: 400 W

Explanation of conformity

Paragraph 4.2 of the EG rules 89/392 CEE

The company IP Dental Division GmbH explains, that the machine, type polish motor

PM 30N

will comply with the EG rules:

IEC 34/EN 60034

As well as the changes CEE 93/68 (1983) as to the electric magnetical compatibility.
Contained safety base rules will answer to and agree with the united standards sizes.

CEE 73/23 (1973)

Instruction

PM 30 N Polishmotor



IP Dent Division GmbH
Pfarrstr. 3
85778 Haimhausen bei München
Deutschland

Telefon:(0049) 08133 - 444951
Telefax:(0049) 08133 - 444953

E-Mail:kontakt@ipdent.de